# **Product information**

# Metal Putty

# PI 24/12/12/2023



### **Description**

Adherent 2-component epoxy resin putty with very good chemical resistance. Ideally suited for permanent quick repairs such as sealing cracks or repairing damaged threaded holes in iron, cast iron and other surfaces. Also suitable for magnesium. After curing, the repaired site can be machined, ground and painted over.

# MONTH TO THE PROPERTY OF THE P

## **Properties**

- rapid curing
- good adhesion
- fills holes, incorrectly drilled holes, cracks and broken off threads
- after curing, suitable for machining
- can be painted over
- high strength

#### Technical data

Base 2-c-epoxy resin metal filled

Color / appearance dark grey

Processing time ca. 3-4 min

Initial strength 8-10 min

Final strength 24 h

Processing temperature +10 to +30 °C
Curing temperature +6 to +30 °C
Combined tension and shear resistance DIN 53283
Mixing ratio 01:01
Shore D hardness 87

Operating temperature range -50°C to +180°C (

short term to

Phone: +49 (0)731/1420-0 Fax: +49 (0)731/1420-82 e-mail: info@liqui-moly.de

+300°C)) ca. 0,05 %

 $\begin{array}{lll} \text{Shrinkage} & \text{ca. 0,05 \%} \\ \text{Thermal conductivity} & \text{0,6 W/m·K} \\ \text{Capacity for electrical} & \text{3 kV/mm} \end{array}$ 

breakdown

Bridging of adhesive gap 15 mm Shelf life in original sealed 36 months

container

Recommended storage 20 °C

temperature

# Areas of application

For sealing cracks in housings, for repairing damaged threaded holes and for securing stud bolts, etc.

# **Application**

To obtain the optimum adhesion, the area to be treated should be free from dirt, grease and oil

residues. For cleaning, we recommend Cleaner and Thinner (part no. 6130). Use a knife to cut off the required length of metal putty from the roll and knead it well.

#### Note:

The two components are mixed together by kneading them together and the metal putty begins to cure. Processing should be performed within 3–4 minutes. The mechanical processing can be performed after approx. 20 minutes.

# Available pack sizes

56 g Blister 6187

D-GB-F-I-E-NL-P

Our information is based on thorough research and may be considered reliable, although not legally binding.